

6356 Manor Lane
South Miami, FL, 33143
Phone 866.593.8444
Fax 305.590.5537
info@tcbiomedix.com



Sterile Isopropyl Alcohol Wipers

Our sterile isopropyl alcohol wipers are 70% saturated and are low linting. Composed of a 55/45 cellulose blend, each 9" x 9" wiper is 65 gm/sq meter, with 30 wipers per zip pack, and 24 packs per case.

- Gamma Irradiated to a 10⁻⁶ Sterility Assurance level
- Saturated with Isopropyl Alcohol, 70% USP
- Alcohol is initially sterilized twice through 0.2 micron Filters, and then terminally sterilized per the Gamma Irradiation highlighted above
- 9" x 9" size
- Extremely low linting, 65gm./sq. meter cellulose blend
- Double bagged for added cleanroom protection



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBWALC30