



IV Administration Set in 92" Length

Our IV Administration Set consists of a universal spike, male swivel luer lock, roller clamp, and drip chamber with 15 micron filter. Each set is individually blister-packaged in quantities of 50 per box and features patient information directly on the exterior of the package.

- 20 drops/ml
- One Y-injection site (6" from distal end)
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15 micron filter located in drip chamber
- 92" length
- DEHP-free
- Latex-free



Product # TCBINF001