

6356 Manor Lane
South Miami, FL, 33143
Phone 866.593.8444
Fax 305.590.5537
info@tcbiomedix.com



IV Administration Set With GVS EasyDrop® Flow Regulator

Our IV Administration Set with Flow Regulator consists of a universal spike, male swivel luer lock, roller clamp, flow regulator, and drip chamber. Each set is individually blister-packaged in quantities of 40 per box and features patient information directly on the exterior of the package.

- 20 drops/ml
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15 micron filter located in drop chamber
- GVS EasyDrop® Flow Regulator
- 92" length
- DEHP-free
- Latex-free



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBINF033G