



## Pharmacy Transfer Tubing (Weighted)

*Our custom Single Lead Pharmacy Fluid Transfer Set features large volume intake tubing, a weighted male luer lock intake connector, and silicone tubing insert.*

- For use with Baxa Repeater® Pump
- Large volume intake tubing
- Weighted male luer lock intake connector
- Silicone tubing Insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Sterile



Product # TCBPTT013

Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.