

6356 Manor Lane, Suite 101  
South Miami, FL, 33143  
Phone 866.593.8444 Fax  
305.590.5537  
info@tcbiomedix.com



Packaged in quantities of  
50 per box

## H Connectors

Our Polypropylene H Transfer connector is designed to connect a male luer or luer slip to syringes, tube sets, repeating syringes, bag ports, and other containers. The connectors are:

- Individually packaged
- Sterile
- Latex-Free
- DEHP-Free

Product # TCBPTTH2

Male Luer Lock/Male Luer Slip to Female Luer Lock

