



Filtered IV Extension Sets

- 12" Tip to Tip
- Priming Volume 4.5 ml
- DEHP-Free and Latex Free
- PES Filter Membrane - 0.2 micron and 1.2 micron
- Slide Clamp
- Standard Bore
- Male and Female Luer Lock
- Individually packaged in box of 50 sets



TCBEXT0.2



TCBEXT1.2

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBEXT0.2
Product # TCBEXT1.2
MKT-024 Rev A