




# TCBINFSP1

## IV Administration Gravity Set

- + Standardize practice from pharmacy admixture to patient infusion.
- + Needle-free admixture and infusion.
- + Reduce potential needlestick injuries to healthcare workers.
- + Allows for a safe and efficient infusion process.



Product #	Length	Drops/mL	Priming Vol.	Features	Quantity
TCBINFSP1	83"	10	22mL (approx.)	<ul style="list-style-type: none"> <li>• Integrated Needle-Free Valve Drip Chamber</li> <li>• Needle-Free Y-site</li> <li>• Swirling Luer with Filter Cap</li> <li>• DEHP-Free and Latex-Free</li> </ul>	50 per box

