



Sterile Lab Gowns for Compounding

Our sterile SMS 40gsm weight white anti-static treated cleanroom gowns are designed for non-hazardous drug compounding. The gowns are low linting and latex-free. They feature a stand-up collar, elastic cuffs and thumb loops. They come individually packaged 1 pc/sterile bag with 50 gowns to a case.

- Low linting
- Thumb loops
- Stand up collar
- ISO-certified breathable fabric
- Latex-Free
- Available in sizes Small, Regular, XL and XXXL



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBA40ST-S
Product # TCBA40ST-R
Product # TCBA40ST-XL
Product # TCBA40ST-XXXL