



## IV Administration Set in 92" Length

*Our IV Administration Set consists of a universal spike, male swivel luer lock, roller clamp, and drip chamber with 15µm filter. Each set is individually blister-packaged in quantities of 50 per box.*

- 20 drops/mL
- One Y-injection site (6" from distal end)
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15µm filter located in drip chamber
- 92" length
- DEHP-Free
- Latex-Free



Product # TCBINF001

Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.