

6356 Manor Lane
South Miami, FL, 33143
Phone 866.593.8444
Fax 305.590.5537
info@tcbiomedix.com



Secondary IV Administration Set

- 42" Length
- Universal Vented/Non-Vented Spike
- 20 drops/ml
- Roller Clamp
- Spin Connector with Hanger
- Priming volume ~11 ml
- DEHP-Free PVC/Latex Free



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBINF042
MKT-023 Rev A