



## IV Administration Set

### With GVS SpeedFlow<sup>®</sup> 0.2 $\mu$ m Filter and GVS EasyDrop<sup>®</sup> Flow Regulator

*Our IV Administration Set with 0.2 $\mu$ m filter and flow regulator consists of a universal spike, male swivel luer lock, roller clamp, drip chamber with 15 $\mu$ m filter, and 0.2 $\mu$ m in-line filter. Each set is individually blister-packaged in quantities of 40 per box.*

- 20 drops/mL
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15 $\mu$ m filter located in drip chamber
- GVS Speedflow<sup>®</sup> 0.2 $\mu$ m Filter
- GVS EasyDrop<sup>®</sup> Flow Regulator
- 92" length
- DEHP-Free
- Latex-Free



Product # TCBINF044G

Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.