



Anesthesia IV Administration Set

Our IV Administration Set consists of a universal spike, male swivel luer lock, roller clamp, and 15 drop/ml drip chamber . Each set is individually blister-packaged in quantities of 50 per box.

- 120" Length with One Needle Free Y-Site, Back Check Valve and Stopcock
- Priming Volume: 19ml (approx.)
- Universal Vented/Non Vented 15 Drop/ml. Spike
- Non-DEHP PVC Fluid Path
- Pyrogen Free
- Latex-Free
- Male Luer Lock at Distal End
- Sterile



Product # TCBINF120SC

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.