



Sterile Isopropyl Alcohol Wipers

Our sterile isopropyl alcohol wipers are 70% saturated and are low linting. Composed of a 55/45 cellulose blend, each wiper is 20cm x 23cm and 65gsm, with 30 wipers per zip pack, and 24 packs per case.

- Gamma irradiated to a 10^{-6} sterility assurance level
- Saturated with isopropyl alcohol, 70% USP
- Alcohol is filtered twice through $0.2\mu\text{m}$ filters and then terminally sterilized
- 20cm x 23cm size
- Extremely low linting, 65gsm cellulose blend
- Double bagged for added cleanroom protection



Product # TCBWALC30

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.