



## Sterile Cleanroom Wipers

*Our sterile cleanroom wipers are sterile gamma irradiated to a 10<sup>-6</sup> sterility assurance level. They are constructed of an extremely low linting 68gsm 55/45 cellulose blend. They are appropriate for critical areas up to and including ISO 5 cleanroom environments. The Cleanroom Wipers feature maximum absorbencies of 20-25 mL, and are available in 9" x 9" and 12" x 12" sizes, and available in different quantities.*

### Product # TCBWIP09

300 wipers per pack, 12 packs per case, a total of 3600 wipers per case

### Product # TCBWIP09SP-20

20 wipers per bag, 15 bags per bundle, 12 bundles per case, a total of 3600 wipers per case

### Product # TCBWIP12

150 wipers per pack, 8 packs per case, a total of 1200 wipers per case

### Product # TCBWIP12SP-20

20 wipers per bag, 8 bags per bundle, 8 bundles per case, a total of 1280 wipers per case



Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBWIP09  
Product # TCBWIP09SP-20  
Product # TCBWIP12  
Product # TCBWIP12SP-20